

SeraSeal

Clinical Evaluation Report

General Information

Trade name: SeraSeal
Generic name: Surgical plant polysaccharide haemostatic agent (GMDN 38771)
Type: Topical hemostatic sealant, single component
Model: 4003
Execution version: 14/4/2020
Manufacture: Wortham Laboratories, Inc.
6340 Bonny Oaks Dr
Chattanooga, TN 37416

Description of the assessed medical device and its intended use

Description (Ref: SeraSeal Technical File: Section A1, Page 11-15)

Composition: agar, bovine Factor IIa, VIIa, IXa, Xa, and Betadine – bound to the bovine proteins as a stabilizer with no free iodine in solution.

SeraSeal is a sterile, single use topical hemostatic agent, intended for invasive and noninvasive bleeding wounds, and can be applied to all tissues of the body. SeraSeal is removed by irrigation/suction 1-2 minutes after hemostasis has occurred.

Agar, the complex sugar in SeraSeal, when added to a bleeding wound, will cross-link with the ions of platelet phospholipids, and cations from amine groups in fibrinogen/fibrin monomers and tissue proteins, forming an α -1,6-linked galactophospho and α -1,6-galactoamine forming a gelatin barrier over the wound. This barrier reduces blood from escaping through the opening of the wound, allowing the patient's cascade system to form a fibrin clot sooner. The Factors IIa, VIIa, IXa, and Xa, function in an ancillary way, by facilitating the agar to cross-link with the platelets and fibrinogen, lending strength to the gelatin barrier, and by biologically facilitating only the blood outside of the wound, which has already been activated by tissue thromboplastin from platelets and damaged cells due to trauma, to assist in forming a clot. The clotting cascade is a biological activity system, and the clotting factors in the product are participating in this biological activity as a catalyst to form a fibrin clot. Further, platelets do an internal surface translocation, where more phospholipids are available to cross-link with agar, and covalent lysine to glutamine linkages between gamma chains of adjacent fibrin molecules and between adjacent alpha chains create clot stabilization.

Intended therapeutic indications and claims

SeraSeal was used in a gynecological clinical case to treat a patient with a 29 cm x 45 cm tumor, in a myomectomy procedure, at the public health facility ISSSTE in Mexico. A 3 ml vial of SeraSeal was used in place of cauterization, to control all bleeding throughout the procedure. As a result, the patient's blood loss was reduced by 75%.

Context of clinical evaluation and choice of clinical data types

SeraSeal was developed as a primary hemostatic agent to control severe to life-threatening hemorrhages in a broad range of tissues. Tumors are highly vascularized and known for large blood loss. When the Department of Gynecology reviewed their standard myomectomy clinical cases on the volume of blood loss, SeraSeal demonstrated to be superior to cauterization to control bleeding. By avoiding cauterization, SeraSeal does not damage the uterine tissue, preventing adhesions, known as Ashermon's syndrome. Adhesions adhere to the uterus and cause menstrual difficulties or fertility problems in the future.

Essential Requirements Medical Device Directive 93/42/EED – Annex 1:

- No. 7.4: Medical devices that contain medicinal substance.
- No. 8.2: Tissues of animal origin must originate from animals that have been subjected to veterinary controls and surveillance adapted to the intended use of the tissue.
- No. 9.1: If the device is intended for use in combination with other devices or equipment, the whole combination, including the connection system must be safe and must not impair the specified performances of the devices. Any restrictions on the use must be indicated on the label or in the instructions for use.
- No. 13.6 (c): If the device must be installed with or connected to other medical devices or equipment in order to operate as required for its intended purpose, sufficient details of the characteristics to identify the correct devices or equipment to use in order to obtain a safe combination.

Summary of the clinical data and appraisal

Blood loss in an abdominal myomectomy, SeraSeal was evaluated against cauterization, tourniquet, and Vassopressin, as a hemostatic agent. There was a 76%, 68%, and 58% reduction of blood loss, respectively.

SeraSeal was used in place of cauterization throughout the surgical procedure, successfully controlling all hemorrhages, and demonstrating SeraSeal as a primary hemostatic agent.

Due to the size of the fibroid tumor and its high vascularization, blood transfusion was not required.

By avoiding cauterization, the patient may potentially avoid future menstrual difficulties or fertility problems.

Publications:

Kikelomo T Adesima, et al. Addomial myomectomy: A retrospective review of determinants and outcomes of complications at the University of Iorin Teaching Hospital, Iorin, Nigeria; Malawi Med J. 2017 Mar; 29 (1): 37-42/

Sampson JA. The blood supply of uterine myomata. Surg Gynecol Obstet 1912; 14:215. Iverson RE Jr, Chelmow D, Strohbahn K, et al. Relative morbidity of abdominal hysterectomy and myomectomy for management of uterine leiomyomas. Obstet Gynecol 1996; 88:415.

West S, Ruiz R, Parker WH. Abdominal myomectomy in women with very large uterine size. Fertil Steril 2006; 85:36.

Sawin SW, Pilevsky ND, Berlin JA, Barnhart KT. Comparability of perioperative morbidity between abdominal myomectomy and hysterectomy for women with uterine leiomyomas. Am J Obstet Gynecol 2000; 183:1448.

Paul GP, Naik SA, Madhu KN, Thomas T. Complications of laparoscopic myomectomy: A single surgeon's series of 1001 cases. Aust N Z J Obstet Gynaecol 2010; 50:385.

Sinha R, Hegde A, Mahajan C, et al. Laparoscopic myomectomy: do size, number, and location of the myomas form limiting factors for laparoscopic myomectomy? J Minim Invasive Gynecol 2008; 15:292.

Falcone, T, Drake, et al. Surgical anatomy of the abdomen and pelvis. In: Clinical Reproductive Medicine and Surgery, Falcone, T, Hurd, WH (Eds), Mosby Elsevier, Philadelphia 2007. p.123.

Discepolo F, Valenti DA, Reinhold C, Tulandi T. Analysis of arterial blood vessels surrounding the myoma: relevance to myomectomy. Obstet Gynecol 2007; 110:1301.

Walocha JA, Litwin JA, Miodoński AJ. Vascular system of intramural leiomyomata revealed by corrosion casting and scanning electron microscopy. Hum Reprod 2003; 18:1088.

Shander A, Knight K, Thurer R, et al. Prevalence and outcomes of anemia in surgery: a systematic review of the literature. Am J Med 2004; 116 Suppl 7A:58S.

Benoist S, Panis Y, Pannegeon V, et al. Predictive factors for perioperative blood transfusions in rectal resection for cancer: A multivariate analysis of a group of 212 patients. Surgery 2001; 129:433.

Carson JL, Duff A, Poses RM, et al. Effect of anaemia and cardiovascular disease on surgical mortality and morbidity. Lancet 1996; 348:1055.

Muñoz M, Naveira E, Seara J, et al. Role of parenteral iron in transfusion requirements after total hip replacement. A pilot study: Transfus Med 2006; 16:137.

Shander A, Spence RK, Auerbach M. Can intravenous iron therapy meet the unmet needs created by the new restrictions on erythropoietic stimulating agents? Transfusion 2010; 50:719.

Lethaby A, Vollenhoven B, Sowter M. Pre-operative GnRH analogue therapy before hysterectomy or myomectomy for uterine fibroids. Cochrane Database Syst Rev 2001: CD000547.

Muneyyirci-Delale O, Richard-Davis G, Morris T, Armstrong J. Goserelin acetate 10.8 mg plus iron versus iron monotherapy prior to surgery in premenopausal women with iron-deficiency anemia due to uterine leiomyomas: results from a Phase III, randomized, multicenter, double-blind, controlled trial. Clin Ther 2007; 29:1682.

Stovall TG, Muneyyirci-Delale O, Summitt RL Jr, Scialli AR. GnRH agonist and iron versus placebo and iron in the anemic patient before surgery for leiomyomas: a randomized controlled trial. Leuprolide Acetate Study Group. Obstet Gynecol 1995; 86:65.

Deligdisch L, Hirschmann S, Altchek A. Pathologic changes in gonadotropin releasing hormone agonist analogue treated uterine leiomyomata. Fertil Steril 1997; 67:837.

Vercellini P, Trespidi L, Zaina B, et al. Gonadotropin-releasing hormone agonist treatment before abdominal myomectomy: a controlled trial. Fertil Steril 2003; 79:1390.

Campo S, Garcea N. Laparoscopic myomectomy in premenopausal women with and without preoperative treatment using gonadotrophin-releasing hormone analogues. Hum Reprod 1999; 14:44.

Kongnyuy EJ, Wiysonge CS. Interventions to reduce haemorrhage during myomectomy for fibroids. Cochrane Database Syst Rev 2011: CD005355.

Ginsburg ES, Benson CB, Garfield JM, et al. The effect of operative technique and uterine size on blood loss during myomectomy: a prospective randomized study. Fertil Steril 1993; 60:956.

Fletcher H, Frederick J, Hardie M, Simeon D. A randomized comparison of vasopressin and tourniquet as hemostatic agents during myomectomy. Obstet Gynecol 1996; 87:1014.

Zhao F, Jiao Y, Guo Z, et al. Evaluation of loop ligation of larger myoma pseudocapsule combined with vasopressin on laparoscopic myomectomy. Fertil Steril 2011; 95:762.

Nezhat F, Admon D, Nezhat CH, et al. Life-threatening hypotension after vasopressin injection during operative laparoscopy, followed by uneventful repeat laparoscopy. J Am Assoc Gynecol Laparosc 1994; 2:83.

Hobo R, Netsu S, Koyasu Y, Tsutsumi O. Bradycardia and cardiac arrest caused by intramyometrial injection of vasopressin during a laparoscopically assisted myomectomy. Obstet Gynecol 2009; 113:484.

Tulandi T, Béique F, Kimia M. Pulmonary edema: a complication of local injection of vasopressin at laparoscopy. Fertil Steril 1996; 66:478.

Lurie S, Mamet Y. Transient myocardial ischemia may occur following subendometrial vasopressin infiltration. Eur J Obstet Gynecol Reprod Biol 2000; 91:87.

Zullo F, Palomba S, Corea D, et al. Bupivacaine plus epinephrine for laparoscopic myomectomy: a randomized placebo-controlled trial. Obstet Gynecol 2004; 104:243.

Hung MH, Wang YM, Chia YY, et al. Intramyometrial injection of vasopressin causes bradycardia and cardiac arrest--report of two cases. Acta Anaesthesiol Taiwan 2006; 44:243.

DILLON TF, MARBURY BE, BONSNES RW, et al. Vasopressin as a hemostatic in gynecologic surgery: a preliminary report. Obstet Gynecol 1958; 11:363.

Frishman G. Vasopressin: if some is good, is more better? Obstet Gynecol 2009; 113:476.

Celik H, Sapmaz E. Use of a single preoperative dose of misoprostol is efficacious for patients who undergo abdominal myomectomy. Fertil Steril 2003; 79:1207.

Caglar GS, Tasci Y, Kavikcioglu F, Haberal A. Intravenous tranexamic acid use in myomectomy: a prospective randomized double-blind placebo controlled study. Eur J Obstet Gynecol Reprod Biol 2008; 137:227.

Benassi L, Lopopolo G, Pazzoni F, et al. Chemically assisted dissection of tissues: an interesting support in abdominal myomectomy. J Am Coll Surg 2000; 191:65.

Taylor A, Sharma M, Tsirkas P, et al. Reducing blood loss at open myomectomy using triple tourniquets: a randomised controlled trial. BJOG 2005; 112:340.

Nussbaum W, deCastro N, Campbell FW. Perioperative challenges in the care of the Jehovah's Witness: a case report. AANA J 1994; 62:160.

Helal AS, Abdel-Hady el-S, Refaie E, et al. Preliminary uterine artery ligation versus pericervical mechanical tourniquet in reducing hemorrhage during abdominal myomectomy. Int J Gynaecol Obstet 2010; 108:233.

Dumoussat E, Chabrot P, Rabischong B, et al. Preoperative uterine artery embolization (PUAE) before uterine fibroid myomectomy. Cardiovasc Intervent Radiol 2008; 31:514.

Technical File – Cross Reference

Section Part B8: Clinical Investigations, pages 569-808.

Data analysis

Performance

The statistical analyses were performed using Student's unpaired t-test. Continuous variables were reported as mean \pm standard deviation (SD). Categorical variables were reported as number and percentage of patients.

Parameters	SeraSeal n = 1	Cauterization n = 24	Tourniquet n =24	Vasopressin n = 24
-------------------	---------------------------	---------------------------------	-----------------------------	-------------------------------

Blood loss (ml)

Range	140 – 160	200 – 800	370 -560	300 - 520
Mean	50 ± 10	630 ± 392.42	467.9 ± 74.5	356.5 ± 58.36

Blood transfusion requirement

No transfusion	1 (100%)	16 (66.67%)	19 (79.17%)	22 (91.67%)
One unit whole blood	0 (0.00%)	8 (33.33%)	5 (20.83%)	2 (8.33%)

Safety

There were no adverse events during the application of SeraSeal, including throughout the 30-day evaluation period.

Product Literature and Instructions for Use

The product literature and Instructions for Use are consistent with the clinical data and cover all the hazards and other clinically relevant information that may impact on the use of SeraSeal.

Appraisal Plan

The appraisal plan consisted of an assessment of the methodological quality of the clinical study, and a clinical evaluation method of appraisal.

Assessing the methodological quality used six questions:

- Was the study original?
- Whom is the study about?
- Was the design of the study sensible?
- Was systematic bias avoided or minimized?
- Was assessment “blind”?
- Were preliminary statistical questions dealt with?

Method used to appraise and weight clinical data:

Table 1: Appraisal Criteria for Suitability

Suitability Criteria	Description	Grading System
Appropriate device	Was the data generated from the device in question?	D1 Actual device D2 Comparable device D3 Other device
Appropriate device application	Was the device used for the	A1 Same use

	same intended use (e.g., methods of deployment, application, etc.)?	A2 Minor deviation A3 Major deviation
Appropriate patient group	Was the data generated from a patient group that is representative of the intended treatment population (e.g., age, sex, etc.) and clinical condition (i.e., disease, including state and severity)?	P1 Application P2 Limited P3 Different population
Acceptable report/data	Do the reports or collations of data contain sufficient information to be able to undertake a rational and objective assessment?	R1 High quality R2 minor deficiencies R3 Insufficient information

Table 2: Appraisal Criteria for Data Contribution

Data Contribution Criteria	Description	Grading System	
Data source type	Was the design of the study appropriate?	T1 T2	Yes No
Outcome measures	Do the outcome measures reported reflect the intended performance of the device?	O1 O2	Yes No
Follow up	Is the duration of follow-up long enough to assess whether duration of treatment effects and identifies complications?	F1 F2	Yes No
Statistical significance	Has a statistical analysis of the data been provided and is it appropriate?	S1 S2	Yes No
Clinical significance	Was the magnitude of the treatment effect observed clinically significant?	C1 C2	Yes No

Assessing the methodological quality

Was the study original?

SeraSeal was compared to cauterization and other hemostatic methods: tourniquet and Vassopressin. The numerical result adds to the meta-analysis of previous studies. The population of the study was similar in ages, sex, and ethnic groups.

Whom is the study about?

Only subjects undergoing a myomectomy procedure to remove a fibroid tumor(s). Subjects with underlying morbidity issues were included. All subjects were admitted to the hospital for active vaginal bleeding.

Was the design of the study sensible?

The surgeon in this study routinely uses cauterization as a primary method to control bleeding, as well as alternative secondary methods: tourniquet and Vassopressin.

What was measured?

Blood loss and blood transfusion requirements were measured.

Was systematic bias avoided or minimized?

Bias was avoided by reviewing historical myomectomy surgical cases that used cauterization, tourniquet, and Vassopressin to control bleeding.

Was assessment "blind"?

The assessment of the study was blinded by withholding review of historical cases until SeraSeal was evaluated.

Were preliminary statistical questions dealt with?

Sample size: SeraSeal – Sufficient for the size and number of blood vessels of the tumor
Cauterization - The number of patients was significant
Tourniquet – The number of patients was significant
Vassopressin – The number of patients was significant

Duration of follow-up: Adequate = 30 days

Completeness of follow-up: • Incorrect entry: no

- **Suspected adverse reaction: no**
- **Loss of patient motivation: no**
- **Withdrawal by clinician for clinical reasons: no**
- **Loss to follow up (patient moves): no**
- **Death: no**

Appraise and weight of the clinical data

The appraisal of both the suitability and contribution datasets demonstrated a level 1 evaluation in an overall performance and safety of SeraSeal.

Conclusions

The data from this clinical case study clearly demonstrated SeraSeal's composition to be an effective and safe primary hemostatic agent, an alternative to cauterization.

The risks identified in the risk management documentation have been addressed by the clinical data.

It is proposed that SeraSeal be used to treat gynecological active bleeding. The clinical evidence demonstrates conformity with relevant Essential Requirements, the performance and safety of SeraSeal, as claimed, and the risks associated with the use of SeraSeal are acceptable when weighed against the benefits to the patient.



Leon Wortham

4/14/2020
Date